## Amendments to the Specification:

Please insert the following table at line 25, page 23 of the specification:

TABLE 1

RNA fragment No	nucleotides	Kd <sub>app</sub> (in nM)	
1, <u>SEQ. ID NO.15</u>	00000000	not bound	
2, <u>SEQ. ID NO.16</u>	បស្តិកក្រុកក្រុកក្រុកក្រុកក្រុកក្រុកក្រុកក្រ	0.97 +/~ 0.19	
3, <u>SEQ. ID NO.17</u>	(AUUU)₂A	not bound	
4a, <u>SEQ. ID NO.18</u>	(AUUU)₃A	1.40 +/- 0.39	
4a, <u>SEQ. ID NO.19</u>	AUUUAUUUAUUUA	not bound	
4b, SEQ. ID NO. 20	AUU <b>UAUUUA</b>	0.77 +/- 0.25	
4c, <u>SEQ. ID NO.21</u>	AUUUAUUUAUUUA	not bound	
4d, <u>SEQ, ID NO.22</u>	AUUU <b>AUUUA</b> UUUA	not bound	
prel.consensus	N N U/C U N N U/C U U/C		
7a, <u>SEQ. ID NO.23</u>	UAAUUUUU	not bound	
7b, <u>SEQ. ID NO.24</u>	UAUAUUUU	not bound	
7c, <u>SEQ. ID NO.25</u>	UAUUUU <u>A</u> UU	not bound	
7d, <u>SEQ. ID NO.26</u>	UAUUUUAU	not bound	
8a, <u>SEQ. ID NO.27</u>	UACUUUUU	not bound	
8b, <u>SEQ, ID NO.28</u>	<u> </u>	not bound	
8c, <u>SEQ, ID NO.29</u>	UAUUUU <b>C</b> UU	not bound	
5, <u>SEQ. ID NO.30</u>	UAUUAUUUU	1.14 +/- 0.24	
6, <u>SEQ. ID NO.31</u>	<u>A</u> AUUUAUUU	1.01 +/- 0.27	
MOTIF	NNUUNNUUU		

Line 7, page 25:

TABLE 1:

57 CD CAC ( )				
Name		Sequence	Position	
Openers SEQ.ID NO. 32	O <sub>1</sub>	AATATAAAATTTAAATATTT	804-823	
SEQ. ID NO. 33	Oz	TAGAGCCCCTAGGGCTTACA	909-928	
Negative Controls SEQ. ID NO. 34	N <sub>1</sub>	AGTGGGAAGCACTTAATTAC	757-775	
SEQ. ID NO. 35	$N_2$	CATAATAATAATATTTTGG	969-950	

Line 10, Page 28:

## f) Quantitative real-time RT PCR

ŔNA is reverse transcribed to cDNA using the TaqMan RT PCR reagents (Applied Biosystems) and random hexamers for priming following standard protocols. Quantitative RT-PCR is performed with SYBR Green detection on an ABI7700 instrument (Applied Biosystems) with IL-2 specific primers SEQ. ID NO. 36: (forward: 5'-TCACCAGGATGCTCACATTTAAGTT-3; reverse: 5'-GGAGTTTGAGTTCTTCTTCTAGAC ACTGA-3'; primers are a gift from F. Kalthoff, Novartis Institute for Biomedical Research Vienna). EF-1 alpha is used as endogenous control for normalization (primers: forward SEQ. ID NO. 37: 5'-TTTGAGACCAGCAAGTACTATGTGACT-3', reverse 5'-TCAGCCTGAGATGTCCCTGTAA -3'). The ΔΔCt method is used for relative quantification of IL-2 mRNA levels (as described by Applied Biosystems). All presented data are averages from at least 5 identical independent samples and representative of at least two independent experiments using cells from different donors.

Applicant submits that no new material is being added and the amended tables and line amendments bring the specification into conformance with the sequence listing submitted in Response to the Notice of Missing Requirements mailed June 10, 2008.

Respectfully submitted,

Aftorney for Applicants

John T. Prince

Reg. No. 43,019

Novartis Institutes for BioMedical Research, Inc. 400 Technology Square Cambridge, MA 02139

(617) 871-3346

Date: Agrast B, 2008

-4-